

NUVASIVE® COHERE® ALIF SYSTEM
INTERVERTEBRAL BODY FUSION DEVICE
INSTRUCTIONS FOR USE

For Symbols Glossary, please refer to

<https://www.nuvasive.com/eifu/symbols-glossary>

DESCRIPTION

The NuVasive® Cohere® ALIF System Intervertebral Body Fusion Device interbody implants are manufactured from PEEK Scoria (Polyether-ether-ketone) and Ti-6Al-4V conforming to ASTM F1472/ISO 5832-3. The bolts (bone screws) are manufactured from Ti-6Al-4V ELI conforming to ASTM F136 and ISO 5832-3. The implants are available in a variety of different shapes and sizes to suit the individual pathology and anatomical conditions of the patient. The Cohere® ALIF System 10° - 20° lordotic cages may be used as a standalone system. The Cohere® ALIF System 25° - 30° lordotic cages may be used with supplemental internal spinal fixation systems (e.g., posterior pedicle screw and rod system) that are cleared by the FDA for use in the lumbar spine.

INDICATIONS FOR USE

The NuVasive® Cohere® ALIF System Intervertebral Body Fusion Device is indicated for spinal fusion procedures in skeletally mature patients. The Cohere® ALIF System Intervertebral Body Fusion Device 10°-20° lordotic cages may be used as a standalone system. The Cohere® ALIF System Intervertebral Body Fusion Device 25°-30° lordotic cages must be used with supplemental internal spinal fixation systems (e.g., posterior pedicle screw and rod system) that are cleared by the FDA for use in the lumbar spine. The System is designed for use with autogenous bone graft, allogeneic bone graft comprised of cancellous and/or corticocancellous bone graft, or a bone void filler as cleared by FDA for use in intervertebral body fusion to facilitate fusion. The devices are to be used in patients who have had at least six months of non-operative treatment.

The NuVasive® Cohere® ALIF System Intervertebral Body Fusion Device is intended for use in interbody fusions in the lumbar spine from L2 to S1, following discectomy in the treatment of symptomatic degenerative disc disease (DDD) or degenerative spondylolisthesis, and/or spinal stenosis at one or two adjacent levels. DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies. The Cohere® ALIF System Intervertebral Body Fusion Device implants can also be used as an adjunct to fusion in patients diagnosed with multilevel degenerative scoliosis and sagittal deformity; however, when used in these patients at multiple levels and for patients with degenerative spondylolisthesis or spinal stenosis at one or two adjacent levels, the Cohere® ALIF System Intervertebral Body Fusion Device must be used with a supplemental internal spinal fixation system (e.g., pedicle screw system) cleared by the FDA for use in the lumbar spine in addition to the integrated screws.

CONTRAINDICATIONS

Contraindications include but are not limited to:

1. Infection, local to the operative site.
2. Signs of local inflammation.
3. Patients with known sensitivity to the materials implanted.
4. Patients who are unwilling to restrict activities or follow medical advice.
5. Patients with inadequate bone stock or quality.
6. Patients with physical or medical conditions that would prohibit beneficial surgical outcome.
7. Prior fusion at the level(s) to be treated.
8. Use with components of other systems.
9. Reuse or multiple uses.

POTENTIAL ADVERSE EVENTS AND COMPLICATIONS

As with any major surgical procedures, there are risks involved in spinal/orthopedic surgery. Infrequent operative and postoperative complications that may result in the need for additional surgeries include: early or late infection; damage to blood vessels, spinal cord or peripheral nerves, epidural hematoma; pulmonary emboli; loss of sensory and/or motor function; pleural effusions, hemothorax, chylothorax, pneumothorax, subcutaneous emphysema, need for chest tube insertion, intercostal neuralgia, rib fracture, diaphragm injury; atelectasis; impotence; permanent pain and/or deformity. Rarely, some complications may be fatal. The treatment of multilevel degenerative scoliosis may be associated with a lower interbody fusion rate compared to one- and two-level interbody fusions.

Potential risks identified with the use of this system, which may require additional surgery, include:

- Bending, fracture or loosening of implant component(s)
- Loss of fixation
- Nonunion or delayed union
- Fracture of the vertebra
- Proximal junctional kyphosis (PJK)
- Neurological, vascular or visceral injury
- Metal sensitivity or allergic reaction to a foreign body
- Tissue reactions including macrophage and foreign body reactions adjacent to implants
- Infection
- Decrease in bone density due to stress shielding
- Degenerative changes or instability of segments adjacent to fused vertebral levels
- Malalignment of anatomical structures (i.e., loss of normal spinal contours or change in height)
- Pain, discomfort or abnormal sensations due to the presence of the device
- Nerve damage due to surgical trauma
- Bursitis
- Dural leak
- Paralysis

WARNINGS, CAUTIONS AND PRECAUTIONS

The subject device is intended for use only as indicated.

The implantation of spinal systems should be performed only by experienced spinal surgeons with specific training in the use of this spinal system because this is a technically demanding procedure presenting a risk of serious injury to the patient.

Correct selection of the implant is extremely important. The potential for success is increased by the selection of the proper size of the implant. While proper selection can minimize risks, the size and shape of human bones present limitations on the size and strength of implants. Metallic and internal fixation devices cannot withstand the activity levels and/or loads equal to those placed on normal, healthy bone. These devices are not designed to withstand the unsupported stress of full weight or load bearing alone.

When implanted at adjacent levels, it is important to select the appropriate length NuVasive® Cohere® ALIF System bolt and confirm trajectory under intraoperative fluoroscopy in order to avoid potential bolt impingement.

Caution must be taken due to potential patient sensitivity to materials. Do not implant in patients with known or suspected sensitivity to the aforementioned materials.

If fewer than 3 bolts are used in the 4-hole construct, or fewer than 2 bolts are used in the 2-hole construct, then the system is intended to be used with additional supplemental fixation (cleared by the FDA) for use in the lumbar spine.

These devices can break when subjected to the increased load associated with delayed union or nonunion. Internal fixation appliances are load-sharing devices that hold bony structures in alignment until healing occurs. If healing is delayed, or does not occur, the implant may eventually loosen, bend, or break. Loads on the device produced by load bearing and by the patient's activity level will dictate the longevity of the implant.

Corrosion of the implant can occur. Implanting metals and alloys in the human body subjects them to a constantly changing environment of salts, acids, and alkalis, which can cause corrosion. Placing dissimilar metals in contact with each other can accelerate the corrosion process, which in turn, can enhance fatigue fractures of implants. Consequently, every effort should be made to use compatible metals and alloys in conjunction with each other.

Patients with previous spinal surgery at the level(s) to be treated may have different clinical outcomes compared to those without a previous surgery.

Based on fatigue testing results, when using the Cohere® ALIF System, the physician should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc., which may impact on the performance of this system.

Notching, striking, and/or scratching of implants with any instrument should be avoided to reduce the risk of breakage.

Additional care should be taken at the lower levels of the lumbar spine due to the obstruction of anatomical structures, such as the iliac crest and iliac vessels, surgical access for the subject device at these levels may not be feasible.

Care should be taken to ensure that all components are ideally fixated prior to closure.

Patient Education: Preoperative instructions to the patient are essential. The patient should be made aware of the limitations of the implant and potential risks of the surgery. The patient should be instructed to limit postoperative activity, as this will reduce the risk of bent, broken or loose implant components. The patient must be made aware that implant components may bend, break or loosen even though restrictions in activity are followed.

Single Use/Do Not Re-Use: Reuse of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure, material degradation, potential leachables, and transmission of infectious agents.

MRI Safety Information: The Cohere® ALIF System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Cohere® ALIF System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Compatibility: Do not use the Cohere® ALIF System with components of other systems. Unless stated otherwise, NuVasive devices are not to be combined with the components of another system.

PREOPERATIVE WARNINGS

1. Only patients that meet the criteria described in the indications should be selected.
2. Patient condition and/or predispositions such as those addressed in the aforementioned contraindications should be avoided.
3. Care should be used in the handling and storage of the Cohere® ALIF implants. The implants should not be scratched or damaged. Implants and instruments should be protected during storage and from corrosive environments.

For Sterile Implants: Assure highly aseptic surgical conditions and use aseptic technique when removing the Cohere® ALIF implant from its packaging. Inspect the implant and packaging for signs of damage, including scratched or damaged devices or damage to the sterile barrier. Do not use the Cohere® ALIF implants if there is any evidence of damage.

4. Refer to Cleaning and Sterilization Instructions below for all non-sterile parts.
5. Care should be used during surgical procedures to prevent damage to the device(s) and injury to the patient.

POST-OPERATIVE WARNINGS

During the postoperative phase it is of particular importance that the physician keeps the patient well informed of all procedures and treatments.

Damage to the weight-bearing structures can give rise to loosening of the components, dislocation and migration as well as to other complications. To ensure the earliest possible detection of such catalysts of device dysfunction, the devices must be checked periodically postoperatively, using appropriate radiographic techniques.

METHOD OF USE

Please refer to the Surgical Technique for this device.

HOW SUPPLIED

The Cohere® ALIF sterile implants are supplied pre-packaged and sterile. The integrity of the packaging should be checked to ensure that the sterility of the contents is not compromised. Remove the device from the packaging using aseptic technique, only after the correct size has been determined. This product should NOT be re-sterilized. Discard after use.

PACKAGING

Packages for both sterile and non-sterile components should be intact upon receipt. Devices should be carefully examined for completeness, and for lack of damage, prior to use. Damaged packages or products should not be used and should be returned to NuVasive.

Instruments provided non-sterile can be single-use or reusable. Discard single-use instruments after use.

Reusable instruments should be reprocessed using instructions provided below.

Implants provided non-sterile are single use and should be sterilized per the instructions provided below.

All implants and instruments provided sterile are intended for single use only. Do not use if package is opened or damaged. This product should not be re-sterilized. Discard single-use instruments after use.

HANDLING OF THE STERILE IMPLANT

- Before removing the implants from the package, make sure that the protective packaging is unopened and undamaged. If the packaging is damaged, the implants have to be considered as **NON-STERILE** and may not be used.
- Upon removal from the package, compare the descriptions on the label with the package contents (product number and size).
- Note the STERILE expiry date. Implants with elapsed STERILE expiry dates have to be considered as **NON-STERILE** and may not be used.
- Take particular care that aseptic integrity is assured during removal of the implant from the inner packaging.
- Open the pouches carefully, beginning from the triangular corner. Take suitable measures to ensure that the implant does not come into contact with objects that could damage its surfaces. Damaged implants must not be used. Use only the recommended instruments for implantation of the implants.

CLEANING AND DECONTAMINATION

All non-sterile instruments must first be thoroughly cleaned using the validated methods prescribed in the NuVasive Cleaning and Sterilization Instructions (doc #9400896) before sterilization and introduction into a sterile surgical field. Contaminated instruments should be wiped clean of visible soil at the point of use, prior to transfer to a central processing unit for cleaning and sterilization. The validated cleaning methods include both manual and automated cleaning. Visually inspect the instruments following performance of the cleaning instructions to confirm there is no visual contamination of the instruments prior to proceeding with sterilization. If possible contamination is present at visual inspection, repeat the cleaning steps. Contaminated instruments should not be used and should be returned to NuVasive. Contact your local representative or NuVasive directly for any additional information related to cleaning of NuVasive surgical instruments.

Instruments with a "D" prefix part number (e.g., DXXXXXXX) may be disassembled. Please refer to the additional disassembly instructions for these instruments.

STERILIZATION

All non-sterile instruments and implants are sterilizable by steam autoclave using standard hospital practices, in addition to NuVasive's validated parameters. In a properly functioning and calibrated steam sterilizer, effective sterilization may be achieved using the parameters prescribed in the NuVasive Cleaning and Sterilization Instructions (doc #9400896).

INFORMATION

To obtain a Surgical Technique Manual or should any information regarding the products or their uses be required, please contact your local representative or NuVasive directly at +1-800-475-9131. You may also email: info@nuvasive.com.

This Instructions for Use document is intended for the US market only.